



<b>DE</b>	Blutdruckmessgerät Gebrauchsanweisung .....	2
<b>EN</b>	Blood pressure monitor Instructions for use .....	17
<b>FR</b>	Tensiomètre Mode d'emploi .....	31
<b>ES</b>	Tensiómetro Manual de instrucciones .....	46
<b>IT</b>	Misuratore di pressione Istruzioni per l'uso .....	61
<b>TR</b>	Bilgisayarlı tansiyon ölçer Kullanım kılavuzu .....	76
<b>RU</b>	Прибор для измерения давления на запястье Инструкция по применению .....	90
<b>PL</b>	Ciśnieniomierz obsługi Instrukcja obsługi .....	106

# ENGLISH

## Contents

1. Getting to know your instrument.....	17	7. Saving, retrieving and deleting results .....	26
2. Important information.....	18	8. Transferring measurements .....	27
3. Unit description .....	22	9. Error messages/trouble-shooting .....	27
4. Prepare measurement .....	23	10. Cleaning and storing the device and cuff.....	28
5. Measuring blood pressure.....	24	11. Specifications.....	28
6. Evaluating results .....	25	12. Warranty/service.....	29

## Included in delivery:

Blood pressure monitor with cuff  
2 x 1.5 V LR03 AAA batteries  
USB cable  
Storage box  
Instructions for use

## Dear Customer,

thank you for choosing a product from our range. Our name stands for high-quality, thoroughly tested products for applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage, beauty and air. Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

Best regards, Your Beurer Team

## 1. Getting to know your instrument

Check that the device packaging has not been tampered with and make sure that all contents are present. Before use, ensure that there is no visible damage to the device or accessories and that all packaging material has been removed. If you have any doubts, do not use the device and contact your retailer or the specified Customer Services address.

The wrist blood pressure monitor is used for non-invasive measurement and monitoring of adults' arterial blood pressure. You can use it to measure your blood pressure quickly and easily, storing the results and displaying the progression of readings. A warning is issued for anyone suffering from cardiac arrhythmia.

## 2. Important information



### Signs and symbols

The following symbols are used in these instructions for use, on the packaging and on the type plate for the device and accessories:

	Attention
	Note Note on important information
	Observe the instructions for use
	Application part, type BF
	Direct current
	Disposal in accordance with the Waste Electrical and Electronic Equipment EC Directive – WEEE
	Dispose of packaging in an environmentally friendly manner
	Manufacturer

Storage / Transport 	Permissible storage and transport temperature and humidity
Operating 	Permissible operating temperature and humidity
	Protect from moisture
	Serial number
	The CE labelling certifies that the product complies with the essential requirements of Directive 93/42/EEC on medical devices
	Certification symbol for products that are exported to the Russian Federation and members of the CIS



### Advice on use

- In order to ensure comparable values, always measure your blood pressure at the same time of day.
- Before every measurement, relax for about five minutes.
- If you want to perform several measurements on the same person, wait five minutes between each measurement.

- Do not take a measurement within 30 minutes after eating, drinking, smoking or exercising.
- Repeat the measurement if you are unsure of the measured value.
- The measurements taken by you are for your information only – they are not a substitute for a medical examination! Discuss the measurements with your doctor, and never base any medical decisions on them (e.g. medicines and their administration)!
- Using the blood pressure monitor outside your home environment or whilst on the move (e.g. whilst travelling in a car, ambulance or helicopter, or whilst undertaking physical activity such as playing sport) can influence the measurement accuracy and cause incorrect measurements.
- Do not use the blood pressure monitor on newborns, pregnant women or patients with preeclampsia.
- In the case of restricted circulation on the arm as a result of chronic or acute vascular diseases (including vascular constriction), the accuracy of the wrist measurement is limited. In this case you should avoid using an upper arm blood pressure monitor.
- This device is not intended for use by people (including children) with restricted physical, sensory or mental skills or a lack of experience and/or a lack of knowledge, unless they are supervised by a person who is responsible for their safety or are instructed by such a person in how to use the device. Supervise children around the device to ensure they do not play with it.
- Cardiovascular diseases may lead to incorrect measurements or have a detrimental effect on measurement accuracy. The same also applies to very low blood pressure, diabetes, circulatory disorders and arrhythmias as well as chills or shaking.
- The blood pressure monitor must not be used in connection with a high-frequency surgical unit.
- Only use the unit on people who have the specified wrist measurement for the device.
- Please note that when inflating, the functions of the limb in question may be impaired.
- During the blood pressure measurement, blood circulation must not be stopped for an unnecessarily long time. If the device malfunctions, remove the cuff from the arm.
- Do not allow sustained pressure in the cuff or frequent measurements. The resulting restriction of the blood flow may cause injury.



- Ensure that the cuff is not placed on an arm in which the arteries or veins are undergoing medical treatment, e.g. intravascular access or therapy, or an arteriovenous (AV) shunt.
- Do not use the cuff on people who have undergone a mastectomy.
- Do not place the cuff over wounds as this may cause further injury.
- Place the cuff on your upper arm only. Do not place the cuff on other parts of the body.
- The blood pressure monitor can only be operated with batteries.
- To conserve the batteries, the monitor switches off automatically if no buttons are pressed for one minute.
- The device is only intended for the purpose described in these instructions for use. The manufacturer is not liable for damage resulting from improper or careless use.

### **Storage and Care**

- The blood pressure monitor is made up of precision electronic components. Accuracy of readings and the instrument's service life depend on careful handling.
  - You should protect the device from impact, moisture, dirt, major temperature fluctuations and direct exposure to the sun's rays.
  - Never drop the device.

- Do not use near strong electromagnetic fields, i.e. keep it away from any radio systems and mobile phones.
  - Only ever use the cuffs provided with the monitor or original replacement cuffs. Otherwise erroneous results will be recorded.
- If the instrument is not used for any length of time, we recommend removing the batteries.

### **Notes on handling batteries**

- If your skin or eyes come into contact with battery fluid, flush out the affected areas with water and seek medical assistance.
-  **Choking hazard!** Small children may swallow and choke on batteries. Store the batteries out of the reach of small children.
- Observe the plus (+) and minus (-) polarity signs.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Protect the batteries from excessive heat.
-  **Risk of explosion!** Never throw batteries into a fire.
- Do not charge or short-circuit batteries.
- If the device is not to be used for a long period, take the batteries out of the battery compartment.
- Use identical or equivalent battery types only .
- Always replace all batteries at the same time.

- Do not use rechargeable batteries.
- Do not disassemble, split or crush the batteries.

### **Repair and disposal**

- Batteries do not belong in domestic refuse. Used batteries should be disposed of at the collection points provided.
- Never open the instrument. If these instructions are not heeded, the warranty will be null and void.
- Never attempt to repair the instrument or adjust it yourself. We can no longer guarantee perfect functioning if you do.
- Repairs may only be performed by Beurer Customer Service or authorized dealers. However, always check the batteries and replace them if necessary prior to making any complaint.
- For environmental reasons, do not dispose of the device in the household waste at the end of its useful life. Dispose of the unit at a suitable local collection or recycling point. Dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.



### **Notes on electromagnetic compatibility**

- The device is suitable for use in all environments listed in these instructions for use, including domestic environments.
- The use of the device may be limited in the presence of electromagnetic disturbances. This could result in issues such as error messages or the failure of the display/ device.
- Avoid using this device directly next to other devices or stacked on top of other devices, as this could lead to faulty operation. If, however, it is necessary to use the device in the manner stated, this device as well as the other devices must be monitored to ensure they are working properly.
- The use of accessories other than those specified or provided by the manufacturer of this device can lead to an increase in electromagnetic emissions or a decrease in the device's electromagnetic immunity; this can result in faulty operation.
- Failure to comply with the above can impair the performance of the device.

### 3. Unit description



1. Display
2. Memory button **MEM**
3. START/STOP button
4. USB interface
5. Touch screen activation switch
6. Battery cover


#### Touch screen activation switch

The device features a touch screen display. To prevent the screen from being activated accidentally, keep the touch screen activation switch set to **OFF** when the device is not in use. To operate the device, move the touch screen activation switch to **ON**. When the touch screen display is touched (START/STOP button or **MEM** button), a beep is output.

**Tip:** You can switch off the device at any time by setting the touch screen activation switch to the **OFF** position.

#### Icons in the display:

1. Time and date
2. Risk indicator
3. Systolic pressure
4. Diastolic pressure
5. Measured pulse
6. „Change battery“ icon
7. Pump up, release air (arrow)
8. Memory display () , average value (R), morning (R<sup>M</sup>), evening (P<sup>M</sup>)
9. Cardiac arrhythmia icon



#### PC interface

The Beurer blood pressure monitor also allows you to transfer your measured values to the PC.

To do this you require a USB cable (included in delivery) and the Beurer “Health Manager” PC software.

The software can be downloaded free of charge from [www.beurer.com/service/download](http://www.beurer.com/service/download)

#### System requirements for the Beurer “Health Manager” computer software

##### 1. Supported operating systems:

- Windows XP SP3
- Windows Vista SP1 or later
- Windows 7
- Windows 7 SP1

## 2. Supported architectures:

- x86 (32 bit)
- x64 (64 bit)


## 3. Hardware requirements:


- Recommended: at least Pentium 1 GHz or faster with at least 1 GB RAM
- Free memory on the primary partition of at least:
  - x86 – 600 MB
  - x64 – 1.5 GB
- Graphic resolution from: 1024 x 768 pixels
- USB port 1.0 or later

## 4. Prepare measurement

### Inserting battery

- Open the battery compartment lid.
- Use exclusively brand-name batteries of the type: 2 x 1.5V Micro (Alkaline Type LR 03). Making absolutely sure that you insert them with the correct polarity as marked. Do not use rechargeable batteries.
- Replace the battery cover carefully.

If the battery change  is continuously illuminated, measurement is no longer possible and you must replace all the batteries. Once batteries have been removed from the device, the time must be reset.



The empty, completely flat batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.

The codes below are printed on batteries containing harmful substances:

Pb = Battery contains lead,

Cd = Battery contains cadmium,

Hg = Battery contains mercury.



### Setting date and time

You should set the date and time without fail. Otherwise, you will not be able to save your measured values correctly with a date and time and to access them again later.

**Tip:** If you press and hold the **MEM** button, you can set the values more quickly.

To set the date and time, proceed as follows:


- Move the touch screen activation switch to the **ON** position.
- Press the START/STOP and **MEM** buttons simultaneously, 24h begins to flash. Select 12h or 24h mode using the **MEM** button. Press START/STOP to confirm your selection. The year display will start to flash. Set the year with the **MEM** button and confirm with START/STOP ①.
- Then set the month, day, hour and minutes and confirm each setting with START/STOP ①.
- Pressing the START/STOP button ① again will switch the display off.



## 5. Measuring blood pressure

Ensure the device is at room temperature before measuring. The measurement can be performed on the left or right wrist.

### Positioning cuff



- Expose your wrist. Making sure that the circulation in the arm is not restricted by any clothes etc. that are too tight. Position the cuff on the inside of your wrist.
- Fasten the cuff with the Velcro fastening so that the upper edge of the monitor is positioned approx. 1 cm below the ball of your thumb.
- The cuff has to be fitted tightly around the wrist but should not constrict it.

**Important:** The instrument should only be operated with the original cuff.

Blood pressure may vary between the right and left wrist, which may mean that the measured blood pressure values are different. Always perform the measurement on the same wrist.

If the values between the two wrists are significantly different, please consult your doctor to determine which wrist should be used for the measurement.

### Correct posture





- Rest for approx. 5 minutes before each measurement. Otherwise there may be divergences.
- You can perform the measurement while sitting or lying. To carry out a blood pressure measurement, make sure you are sitting comfortably with your arms and back leaning on something. Do not cross your legs. Place your feet flat on the ground. Make sure to rest your arm and move it. Always make sure that the cuff is at heart level. Otherwise significant deviations can occur. Relax your arm and the palm of your hand.
- In order not to distort the result, it is important to keep still during the measurement and not talk.



### Measuring blood pressure

- Move the touch screen activation switch to the **ON** position.
- Switch on the blood pressure monitor with the START/STOP button ①. After the full-screen display, the most recently used user memory appears (U1 or U2). To change the user memory, press the **MEM** button and confirm your selection with the START/STOP button ①. If no but-


ton is pressed, the most recently used user memory is automatically used after 5 seconds.

- Before the measurement, the last saved test result is briefly displayed. If there is no measurement in the memory, the instrument always displays the value  $\square$ .
- The cuff automatically inflates. Cuff air pressure is released slowly. If a tendency towards high blood pressure is already detectable, the cuff is pumped up again and cuff pressure increased further. As soon as a heart rate is detected, the heart rate symbol  is displayed.
- After the pressure has been completely reduced, the systolic pressure, diastolic pressure and pulse readings are displayed.
- You can cancel the measurement at any time by pressing the START/STOP button  or by pushing the touch screen activation button to **OFF**.
-  appears if it has not been possible to perform the measurement properly. Observe the section in these instructions for use on error messages/troubleshooting and repeat the measurement.
- The test result is saved automatically.
- To switch off the device, press the START/STOP button  or push the touch screen activation switch to **OFF**. If you forget to switch off the device, it switches off automatically after approx. 1 minute.



Wait at least 5 minutes before taking another measurement!

## 6. Evaluating results

### Cardiac arrhythmia:

This instrument can identify possible cardiac arrhythmia disorders during measurement and if necessary indicates the measurement with the flashing icon .

This may be an indicator for arrhythmia. Arrhythmia is a condition where the heart rhythm is abnormal as a result of defects in the bioelectrical system controlling the heart beat. The symptoms (omitted or premature heart beats, slow or excessively fast heart rate) may be caused, among other things, by heart disease, age, physical predisposition, excessive use of stimulants, stress or lack of sleep. Arrhythmia can only be ascertained through examination by your doctor.

Repeat the measurement if the flashing icon , is displayed after the measurement. Please note that you should rest for 5 minutes between measurements and not talk or move during the measurement. If the icon , appears often, please contact your doctor. Any self-diagnosis and treatment based on the test results may be dangerous. It is vital to follow your doctor's instructions.

### Risk indicator:

The measurements can be classified and evaluated in accordance with the following table.

However, these standard values serve only as a general guideline, as the individual blood pressure varies in different people and different age groups etc.

It is important to consult your doctor regularly for advice. Your doctor will tell you your individual values for normal blood pressure as well as the value above which your blood pressure is classified as dangerous.

The classification on the display and the scale on the unit show which category the recorded blood pressure values fall into. If the values of systole and diastole fall into two different categories (e.g. systole in the 'High normal' category and diastole in the 'Normal' category), the graphical classification on the device always shows the higher category; for the example given this would be 'High normal'.

Blood pressure value category	Systole (in mmHg)	Diastole (in mmHg)	Action
Setting 3: severe hypertension	≥ 180	≥ 110	seek medical attention
Setting 2: moderate hypertension	160–179	100–109	seek medical attention
Setting 1: mild hypertension	140–159	90–99	regular monitoring by doctor


Blood pressure value category	Systole (in mmHg)	Diastole (in mmHg)	Action
High normal	130–139	85–89	regular monitoring by doctor
Normal	120–129	80–84	self-monitoring
Optimal	< 120	< 80	self-monitoring

Source: WHO, 1999 (World Health Organization)

## 7. Saving, retrieving and deleting results


- The results of each successful measurement are stored together with date and time. With more than 60 items of measured data, the earliest items of data measured are lost.
- Move the touch screen activation switch to the **ON** position.
- Make a selection with the **MEM** button, and then confirm the desired user memory with the START/STOP button **ⓘ**. If you press the **MEM** button again, the average value  $\bar{P}$  of all the stored measured values in the user memory will be displayed. If you press the **MEM** button again, the average value of the last 7 days for the morning measurement will be displayed (Morning: 5 a.m. – 9 a.m., display  $\bar{P}_m$ ). If you press the **MEM** button again, the average value of the last 7 days for the evening measurement will be displayed (Evening: 6 p.m. – 8 p.m., display

Pf)). If you continue to press the **MEM** button, the most recent individual measured values with date and time are displayed in turn.

- You can delete the memory by pressing and holding the **MEM** button for 3 seconds. All the values in the current user memory are deleted after three beeps are output.
- To switch off the device, press the **MEM** button again or the START/STOP button  or push the touch screen activation switch to **OFF**.
- If you forget to switch off the device, it will switch off automatically after 2 minutes.

## 8. Transferring measurements

Connect the blood pressure monitor to your PC using the USB cable.

 No data transfer may be launched whilst performing a measurement.

**PC** is shown on the display. Begin the data transfer in the “HealthManager” PC software. During the data transfer, an animation is shown on the display. A successful data transfer is displayed as in figure 1.



figure 1

If the data transfer is unsuccessful, an error message appears as in figure 2. In this case, interrupt the PC connection and start the data transfer again.

After 30 seconds of not being in use or if communication with the PC is interrupted, the blood pressure monitor switches itself off automatically.




figure 2

## 9. Error messages/trouble-shooting

In case of faults, the **Er** message appears in the display.

Error messages can occur when

- systolic or diastolic pressure could not be measured (**Er 1** or **Er 2** appears on the display),
- systolic or diastolic pressure was outside the measurement range (**Hi** or **Lo** appears on the display),
- the cuff is fastened too tightly or loosely (**Er 3** or **Er 4** appears on the display),
- the pump pressure is higher than 300 mmHg (**Er 5** appears on the display),
- pumping up takes longer than 160 seconds (**Er 6** appears on the display),
- there is a system or device error (**Er A**, **Er 0**, **Er 7** or **Er 8** appears on the display),
- the batteries are almost empty ,
- the data could not be sent to the PC (**PC Er** appears in the display).

In such cases, repeat the measurement. Ensure that you do not move or speak during the measurement. If necessary, reinsert or replace the batteries.

### **Technical alarm – description**

Should the recorded blood pressure (systolic or diastolic) lie outside the limits specified in the section “Technical specifications”, the technical alarm will appear on the display indicating either “**Hi**” or “**Lo**”. In such cases, you should seek medical assistance and check the accuracy of your procedure.

The limit values for the technical alarm are factory set and cannot be adjusted or deactivated. These alarm limit values are accorded second priority under the standard IEC 60601-1-8.

The technical alarm is a non-locking alarm and must not be reset. The signal shown on the display will disappear automatically after about 8 seconds.


## **10. Cleaning and storing the device and cuff**

- Clean the device and cuff carefully using a slightly damp cloth only.
- Do not use any cleaning agents or solvents.
- Under no circumstances hold the device and cuff under water, as this can cause liquid to enter and damage the device and cuff.

- If you store the device and cuff, do not place heavy objects on the device and cuff. Remove the batteries.

## **11. Specifications**

Model no.	BC 58
Measurement method	Oscillometric, non-invasive blood pressure measurement on the wrist
Measurement range	Cuff pressure 0–300 mmHg, systolic 60–260 mmHg, diastolic 40–199 mmHg, Pulse 40–180 beats/minute
Display accuracy	Systolic $\pm 3$ mmHg, diastolic $\pm 3$ mmHg, pulse $\pm 5\%$ of the value shown
Measurement inaccuracy	Max. permissible standard deviation according to clinical testing: systolic 8 mmHg/diastolic 8 mmHg
Memory	2 x 60 memory spaces
Dimensions	L 90 mm x W 68 mm x H 30 mm
Weight	Approx. 158 g (without batteries, with cuff)
Cuff size	140 to 195 mm
Permissible operating conditions	+5 °C to +40 °C, $\leq 90\%$ relative air humidity (non-condensing)

Permissible storage conditions	-20 °C to +55 °C, ≤95 % relative air humidity, 800–1050 hPa ambient pressure
Power supply	2 x 1,5V  AAA batteries
Battery life	For approx. 300 measurements, depending on the blood pressure level and/or pump pressure
Classification	Internal supply, IPX0, no AP or APG, continuous operation, type BF applied part

The serial number is located on the device or in the battery compartment.

Technical information is subject to change without notification to allow for updates.

- This unit is in line with European Standard EN 60601-1-2 (In accordance with CISPR 11, IEC 61000-4-2, IEC 61000-4-3, IEC 61000-4-6 and IEC 61000-4-8) and is subject to particular precautions with regard to electromagnetic compatibility (EMC). Please note that portable and mobile HF communication systems may interfere with this unit.
- This device is in line with the EU Medical Devices Directive 93/42/EEC, the “Medizinproduktegesetz” (German Medical Devices Act) and the standards EN 1060-1

(non-invasive sphygmomanometers, Part 1: General requirements), EN 1060-3 (non-invasive sphygmomanometers, Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems) and IEC 80601-2-30 (Medical electrical equipment – Part 2–30: Particular requirements for the safety and essential performance of automated non-invasive blood pressure monitors).

- The accuracy of this blood pressure monitor has been carefully checked and developed with regard to a long useful life. If using the device for commercial medical purposes, it must be regularly tested for accuracy by appropriate means. Precise instructions for checking accuracy may be requested from the service address.

## 12. Warranty/service

In case of a claim under the warranty please contact your local dealer or the local representation which is mentioned in the list “service international”.

In case of returning the unit please add a copy of your receipt and a short report of the defect.

The following warranty terms shall apply:

1. The warranty period for BEURER products is either 5 years or- if longer- the country specific warranty period from date of Purchase.

In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.

2. Repairs (complete unit or parts of the unit) do not extend the warranty period.
3. The warranty shall not be valid for damages because of
  - a. improper treatment, e.g. nonobservance of the user instructions.
  - b. repairs or tampering by the customer or unauthorised third parties.
  - c. transport from the manufacturer to the consumer or during transport to the service centre.
  - d. The warranty shall not be valid for accessories which are subject to normal wear and tear (cuff, batteries etc.).
4. Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.